

WHEN ONLY THE BEST WILL DO

Swingle Ortho Professional Kit

SAFE AND PRECISE INTERPROXIMAL STRIPPING

- Fast & efficient opening of contact points
- Controlled, step by step interproximal reduction
- Contouring, finishing and polishing



Kit Includes:

Contra Angle Handpiece (with/without light),

File Ejector, Spray Canal Cleaning Needle

6 Ortho Strips (08µm, 80µm, 60µm, 40µm, 25µm, 15µm)

Interproximal Stripping Protocol



1. Opening of the contact points - Opener 08µm



2. Creating Space and contouring - Extra Coarse 80µm, Coarse 60µm, Medium 40µm



3. Finishing/Polishing Fine 25µm Polishing 15µm

To order call us now on **0117 975 5533**
or visit **www.tocdental.com**

